

For Research Use Only

LEPROT Polyclonal antibody

Catalog Number: 20304-1-AP



Basic Information

Catalog Number: 20304-1-AP	GenBank Accession Number: BC011027	Purification Method: Antigen affinity purification
Size: 650 µg/ml	GeneID (NCBI): 54741	Recommended Dilutions: WB 1:200-1:1000
Source: Rabbit	UNIPROT ID: O15243	
Isotype: IgG	Full Name: leptin receptor overlapping transcript	
Immunogen Catalog Number: AG14126	Calculated MW: 131 aa, 14 kDa	
	Observed MW: 16 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : pig heart tissue,
Species Specificity: human, pig	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



pig heart tissue was subjected to SDS PAGE followed by western blot with 20304-1-AP (LEPROT antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.