

TSPAN32 Polyclonal antibody

Catalog Number: 20320-1-AP

Basic Information

Catalog Number:

20320-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14174

GenBank Accession Number:

BC016693

GeneID (NCBI):

10077

UNIPROT ID:

Q96QS1

Full Name:

tetraspanin 32

Calculated MW:

320 aa, 35 kDa

Observed MW:

33 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse spleen tissue,

Background Information

Storage

Storage:


Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

Selected Validation Data



mouse spleen tissue were subjected to SDS PAGE followed by western blot with 20320-1-AP (TSPAN32 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.