

For Research Use Only

# MPDU1 Polyclonal antibody, PBS Only

Catalog Number:20321-1-PBS



## Basic Information

Catalog Number:

20321-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14178

GenBank Accession Number:

BC001898

GeneID (NCBI):

9526

UNIPROT ID:

O75352

Full Name:

mannose-P-dolichol utilization defect

1

Calculated MW:

247 aa, 27 kDa

Observed MW:

20 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 20321-1-AP (MPDU1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20321-1-PBS in a different storage buffer formulation.