

For Research Use Only

HRH4-Specific Polyclonal antibody

Catalog Number: 20362-1-AP



Basic Information

Catalog Number:

20362-1-AP

Size:

727 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_021624

GeneID (NCBI):

59340

UNIPROT ID:

Q9H3N8

Full Name:

histamine receptor H4

Calculated MW:

45 kDa

Observed MW:

~50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse liver tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 20362-1-AP (HRH4-Specific Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.