

For Research Use Only

# RAPSN Polyclonal antibody

Catalog Number: 20453-1-AP

1 Publications



## Basic Information

Catalog Number:

20453-1-AP

Size:

520 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13896

GenBank Accession Number:

BC004196

GeneID (NCBI):

5913

UNIPROT ID:

Q13702

Full Name:

receptor-associated protein of the synapse

Calculated MW:

412 aa, 46 kDa

Observed MW:

43 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Positive Controls:

WB : U-937 cells, NIH/3T3 cells, Y79 cells

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Yue He	37795532	J Appl Physiol (1985)	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



U-937 cells were subjected to SDS PAGE followed by western blot with 20453-1-AP (RAPSN Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.