

For Research Use Only

TAPBPL Polyclonal antibody

Catalog Number: 20522-1-AP



Basic Information

Catalog Number: 20522-1-AP	GenBank Accession Number: BC015017	Purification Method: Antigen affinity purification
Size: 250 µg/ml	GeneID (NCBI): 55080	
Source: Rabbit	UNIPROT ID: Q9BX59	
Isotype: IgG	Full Name: TAP binding protein-like	
Immunogen Catalog Number: AG14371	Calculated MW: 468 aa, 50 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

TAPBPL, also named as TAPBP-R and TAPASIN-R, is a member of the variable-constant Ig superfamily that links MHC class I molecules to the transporter associated with antigen processing in the endoplasmic reticulum.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data