

For Research Use Only

FAM107B Polyclonal antibody, PBS Only

Catalog Number: 20527-1-PBS



Basic Information

Catalog Number:

20527-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14389

GenBank Accession Number:

BC004872

GeneID (NCBI):

83641

UNIPROT ID:

Q9H098

Full Name:

family with sequence similarity 107,
member B

Calculated MW:

131 aa, 16 kDa

Observed MW:

38-40 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 20527-1-AP (FAM107B antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20527-1-PBS in a different storage buffer formulation.