

For Research Use Only

# HIST1H3D Polyclonal antibody, PBS Only

Catalog Number: 20532-1-PBS



## Basic Information

Catalog Number:

20532-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14394

GenBank Accession Number:

BC031333

GeneID (NCBI):

8351

UNIPROT ID:

P68431

Full Name:

histone cluster 1, H3d

Calculated MW:

136 aa, 15 kDa

Observed MW:

15-25 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 20532-1-AP (HIST1H3D Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20532-1-PBS in a different storage buffer formulation.