

For Research Use Only

IKK Alpha Polyclonal antibody

Catalog Number: 20565-1-AP

3 Publications



Basic Information

Catalog Number:

20565-1-AP

Size:

193 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001278

GeneID (NCBI):

1147

UNIPROT ID:

O15111

Full Name:

conserved helix-loop-helix
ubiquitous kinase

Calculated MW:

85 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human, mouse

Background Information

Notable Publications

| Author | Pubmed ID | Journal | Application |
|----------------|-----------|---------------------|-------------|
| Rui Wang | 33831690 | Phytomedicine | WB |
| X U Baogui | 36083478 | J Tradit Chin Med | WB |
| Xiao-Long Wang | 34517275 | Int Immunopharmacol | WB |

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

