

For Research Use Only

# PI3 Kinase p110 Delta Polyclonal antibody



Catalog Number: 20585-1-AP

## Basic Information

Catalog Number:

20585-1-AP

Size:

220 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_005026

GeneID (NCBI):

5293

UNIPROT ID:

O00329

Full Name:

phosphoinositide-3-kinase, catalytic, delta polypeptide

Calculated MW:

119 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data