

For Research Use Only

ZNF486 Polyclonal antibody

Catalog Number: 20599-1-AP



Basic Information

Catalog Number:

20599-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14567

GenBank Accession Number:

BC008936

GeneID (NCBI):

90649

UNIPROT ID:

Q96H40

Full Name:

zinc finger protein 486

Calculated MW:

463 aa, 54 kDa

Observed MW:

54 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : COLO 320 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



COLO 320 cells were subjected to SDS PAGE followed by western blot with 20599-1-AP (ZNF486 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.