

For Research Use Only

NGF Polyclonal antibody

Catalog Number: 20600-1-AP

2 Publications



Basic Information

Catalog Number:

20600-1-AP

Size:

180 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14569

GenBank Accession Number:

BC032517

GeneID (NCBI):

4803

UNIPROT ID:

P01138

Full Name:

nerve growth factor (beta polypeptide)

Calculated MW:

241 aa, 27 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human

Cited Species:

rat

Background Information

Notable Publications

| Author | Pubmed ID | Journal | Application |
|-------------|-----------|-------------|-------------|
| C Hu | 29377496 | Eur J Pain | WB |
| Zunhua Yang | 36638655 | Bioorg Chem | IHC |

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

