

For Research Use Only

SLC22A15 Polyclonal antibody, PBS Only

Catalog Number: 20626-1-PBS



Basic Information

Catalog Number:

20626-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14547

GenBank Accession Number:

BC026358

GeneID (NCBI):

55356

UNIPROT ID:

Q8IZD6

Full Name:

solute carrier family 22, member 15

Calculated MW:

535aa, 59 kDa; 547aa, 61 kDa

Observed MW:

61-66 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 20626-1-AP (SLC22A15 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20626-1-PBS in a different storage buffer formulation.