

For Research Use Only

ACF1/BAZ1A Polyclonal antibody

Catalog Number: 20691-1-AP



Basic Information

Catalog Number:

20691-1-AP

Size:

230 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_013448

GeneID (NCBI):

11177

UNIPROT ID:

Q9NRL2

Full Name:

bromodomain adjacent to zinc finger domain, 1A

Calculated MW:

179 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

BAZ1A, also named as ACF1, WCRF180 and hWALp1, belongs to the WAL family. It may play a role in transcriptional regulation. BAZ1A is involved in the formation or maintenance of heterochromatin playing a critical role in developmental control. The antibody is specific to BAZ1A.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data