

For Research Use Only

ALG6 Polyclonal antibody

Catalog Number: 20753-1-AP



Basic Information

Catalog Number:

20753-1-AP

Size:

153 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14176

GenBank Accession Number:

BC001253

GeneID (NCBI):

29929

UNIPROT ID:

Q9Y672

Full Name:

asparagine-linked glycosylation 6,
alpha-1,3-glucosyltransferase
homolog (S. cerevisiae)

Calculated MW:

507 aa, 58 kDa

Observed MW:

58 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 20753-1-AP (ALG6 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.