

For Research Use Only

TBC1D7 Polyclonal antibody, PBS Only

Catalog Number:20769-1-PBS



Basic Information

Catalog Number:

20769-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14756

GenBank Accession Number:

BC007054

GeneID (NCBI):

51256

UNIPROT ID:

Q9PON9

Full Name:

TBC1 domain family, member 7

Calculated MW:

293 aa, 34 kDa

Observed MW:

30 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human, mouse

Background Information

TBC1D7 is a member of the TBC-domain containing protein family. The encoded protein functions as a subunit of the tuberous sclerosis TSC1-TSC2 complex which plays a role in the regulation of cellular growth and differentiation. Mutations in this gene have been associated with autosomal recessive megalencephaly.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




human placenta tissue were subjected to SDS PAGE followed by western blot with 20769-1-AP (TBC1D7 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20769-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 20769-1-AP (TBC1D7 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20769-1-PBS in a different storage buffer formulation.



IP result of anti-TBC1D7 (IP:20769-1-AP, 4ug; Detection:20769-1-AP 1:500) with Neuro-2a cells lysate 1440 ug. This data was developed using the same antibody clone with 20769-1-PBS in a different storage buffer formulation.



IP result of anti-TBC1D7 (IP:20769-1-AP, 4ug; Detection:20769-1-AP 1:500) with HeLa cells lysate 1200 ug. This data was developed using the same antibody clone with 20769-1-PBS in a different storage buffer formulation.