

For Research Use Only

TTI2 Polyclonal antibody

Catalog Number: 20779-1-AP



Basic Information

Catalog Number:

20779-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14776

GenBank Accession Number:

BC007387

GeneID (NCBI):

80185

UNIPROT ID:

Q6NXR4

Full Name:

chromosome 8 open reading frame 41

Calculated MW:

508 aa, 57 kDa

Observed MW:

57 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, HEK-293T cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 20779-1-AP (TTI2 Antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.