

For Research Use Only

PDCD2L Polyclonal antibody

Catalog Number: 20781-1-AP



Basic Information

Catalog Number:

20781-1-AP

Size:

820 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14795

GenBank Accession Number:

BC006146

GeneID (NCBI):

84306

UNIPROT ID:

Q9BRP1

Full Name:

programmed cell death 2-like

Calculated MW:

358 aa, 39 kDa

Observed MW:

39 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse colon tissue, A549 cells, HCT 116 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 20781-1-AP (PDCD2L antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



mouse colon tissue were subjected to SDS PAGE followed by western blot with 20781-1-AP (PDCD2L Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.