

For Research Use Only

EVA1B Polyclonal antibody

Catalog Number: 20793-1-AP



Basic Information

Catalog Number:

20793-1-AP

Size:

550 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14578

GenBank Accession Number:

BC006241

GeneID (NCBI):

55194

UNIPROT ID:

Q9NVM1

Full Name:

family with sequence similarity 176, member B

Calculated MW:

165 aa, 18 kDa

Observed MW:

25-26 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : human placenta tissue,

Background Information

EVA1B, also known as FAM176B or C1orf78, is a single-pass membrane protein that belongs to the FAM176 family. EVA1B is an important paralog of the EVA1A and is associated with the tumor immune microenvironment and clinical prognosis

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human placenta tissue was subjected to SDS PAGE followed by western blot with 20793-1-AP (EVA1B antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.