

For Research Use Only

ZNF385D Polyclonal antibody

Catalog Number: 20806-1-AP



Basic Information

Catalog Number:

20806-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14746

GenBank Accession Number:

BC007212

GeneID (NCBI):

79750

UNIPROT ID:

Q9H6B1

Full Name:

zinc finger protein 385D

Calculated MW:

395 aa, 42 kDa

Observed MW:

42 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse kidney tissue, human testis tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse kidney tissue were subjected to SDS PAGE followed by western blot with 20806-1-AP (ZNF385D antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.