

For Research Use Only

# KAZALD1 Polyclonal antibody

Catalog Number: 20815-1-AP



## Basic Information

Catalog Number:

20815-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14796

GenBank Accession Number:

BC007758

GeneID (NCBI):

81621

UNIPROT ID:

Q96182

Full Name:

Kazal-type serine peptidase inhibitor domain 1

Calculated MW:

304 aa, 33 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IHC 1:50-1:500

## Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human

**Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (\*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

Positive Controls:

IHC: human heart tissue, human skeletal muscle tissue, human small intestine tissue

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human heart tissue slide using 20815-1-AP (KAZALD1 Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human heart tissue slide using 20815-1-AP (KAZALD1 Antibody) at dilution of 1:200 (under 40x lens).