

For Research Use Only

SLC47A1 Polyclonal antibody

Catalog Number: 20898-1-AP **3 Publications**



Basic Information

Catalog Number: 20898-1-AP	GenBank Accession Number: BC010661	Purification Method: Antigen affinity purification
Size: 550 µg/ml	GeneID (NCBI): 55244	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	UNIPROT ID: Q96FL8	
Isotype: IgG	Full Name: solute carrier family 47, member 1	
Immunogen Catalog Number: AG15013	Calculated MW: 570 aa, 62 kDa	
	Observed MW: 62 kDa, 32 and 50 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse heart tissue, rat heart tissue
Cited Applications: WB	
Species Specificity: human, mouse, rat	
Cited Species: human, rat, mouse	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xian-Mei Song	33952086	Xenobiotica	WB
Muhammad Erfan Uddin	35955741	Int J Mol Sci	WB
Shao-An Wang	38140768	Mol Oncol	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 20898-1-AP (SLC47A1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.