

For Research Use Only

# BEND5 Polyclonal antibody

Catalog Number: 20931-1-AP

Featured Product

2 Publications



## Basic Information

Catalog Number:

20931-1-AP

Size:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14930

GenBank Accession Number:

BC007932

GeneID (NCBI):

79656

UNIPROT ID:

Q7L4P6

Full Name:

BEN domain containing 5

Calculated MW:

421 aa, 48 kDa

Observed MW:

48 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IP

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

WB : mouse liver tissue, Raji cells

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Yanzhu Shi	35844785	Int J Biol Sci	WB,IP
Yuting Ma	37197545	J Thorac Dis	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 20931-1-AP (BEND5 Antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.