

For Research Use Only

# SAP30BP Polyclonal antibody

Catalog Number: 20943-1-AP **1 Publications**



## Basic Information

<b>Catalog Number:</b> 20943-1-AP	<b>GenBank Accession Number:</b> BC007592	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 450 µg/ml	<b>GeneID (NCBI):</b> 29115	<b>Recommended Dilutions:</b> WB 1:200-1:1000
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9UHR5	
<b>Isotype:</b> IgG	<b>Full Name:</b> SAP30 binding protein	
<b>Immunogen Catalog Number:</b> AG15051	<b>Calculated MW:</b> 308 aa, 34 kDa	
	<b>Observed MW:</b> 35-40 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : mouse spleen tissue,
<b>Cited Applications:</b> WB, IP, IF	
<b>Species Specificity:</b> human, mouse	
<b>Cited Species:</b> human	

## Background Information

SAP30BP, also named as SAP30-binding protein, is a 308 amino acid protein, which belongs to the HCNGP family. SAP30BP localizes in Nucleus and Interacts with SAP30. DDX54 Induces cell death and may act as a transcriptional corepressor of a gene related to cell survival. SAP30BP may be involved in the regulation of beta-2-microglobulin genes.

## Notable Publications

Author	Pubmed ID	Journal	Application
Changshou Wang	38059508	EMBO J	WB, IP, IF

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



mouse spleen tissue were subjected to SDS PAGE followed by western blot with 20943-1-AP (SAP30BP Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.