

For Research Use Only

# SAP30BP Polyclonal antibody

Catalog Number: 20943-1-AP

Featured Product

1 Publications



## Basic Information

Catalog Number:

20943-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15051

GenBank Accession Number:

BC007592

GeneID (NCBI):

29115

UNIPROT ID:

Q9UHR5

Full Name:

SAP30 binding protein

Calculated MW:

308 aa, 34 kDa

Observed MW:

35-40 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IP, IF

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

WB : mouse spleen tissue,

## Background Information

SAP30BP, also named as SAP30-binding protein, is a 308 amino acid protein, which belongs to the HCNGP family. SAP30BP localizes in Nucleus and Interacts with SAP30. DDX54 Induces cell death and may act as a transcriptional corepressor of a gene related to cell survival. SAP30BP may be involved in the regulation of beta-2-microglobulin genes.

## Notable Publications

Author	Pubmed ID	Journal	Application
Changshou Wang	38059508	EMBO J	WB,IP,IF

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse spleen tissue were subjected to SDS PAGE followed by western blot with 20943-1-AP (SAP30BP Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.