

For Research Use Only

# Advillin Polyclonal antibody

Catalog Number: 20956-1-AP

1 Publications



## Basic Information

Catalog Number:

20956-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15114

GenBank Accession Number:

BC111730

GeneID (NCBI):

10677

UNIPROT ID:

O75366

Full Name:

advillin

Calculated MW:

819 aa, 92 kDa

Observed MW:

92 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

IF

Species Specificity:

human, mouse

Cited Species:

mouse

Positive Controls:

WB : mouse skin tissue, mouse small intestine tissue

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Charlotte Decourt	37068088	PLoS Biol	IF

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse skin tissue were subjected to SDS PAGE followed by western blot with 20956-1-AP (Advillin antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.