

For Research Use Only

FAM55A Polyclonal antibody, PBS Only

Catalog Number:20996-1-PBS



Basic Information

Catalog Number:

20996-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15021

GenBank Accession Number:

BC029049

GeneID (NCBI):

120400

UNIPROT ID:

Q8N323

Full Name:

family with sequence similarity 55, member A

Calculated MW:

547 aa, 63 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, Indirect ELISA

Species Specificity:

human

Background Information

FAM55A, also known as NXPE1, belongs to the NXPE family. FAM55A is a poorly characterized gene implicated in cellular processes such as cell cycle regulation and neuronal development.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded colon cancer slide using 20996-1-AP (FAM55A antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 20996-1-PBS in a different storage buffer formulation.