

For Research Use Only

C1orf92 Polyclonal antibody, PBS Only

Catalog Number: 21031-1-PBS



Basic Information

Catalog Number:

21031-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15344

GenBank Accession Number:

BC033790

GeneID (NCBI):

149499

UNIPROT ID:

Q8N4P6

Full Name:

chromosome 1 open reading frame 92

Calculated MW:

559 aa, 62 kDa

Observed MW:

62 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 21031-1-AP (C1orf92 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 21031-1-PBS in a different storage buffer formulation.