

For Research Use Only

ANKMY2 Polyclonal antibody

Catalog Number: 21065-1-AP



Basic Information

Catalog Number:

21065-1-AP

Size:

280 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15185

GenBank Accession Number:

BC035353

GeneID (NCBI):

57037

UNIPROT ID:

Q8IV38

Full Name:

ankyrin repeat and MYND domain containing 2

Calculated MW:

441 aa, 49 kDa

Observed MW:

49 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells, HepG2 cells, mouse liver tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 21065-1-AP (ANKMY2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.