

For Research Use Only

# ADORA2B Polyclonal antibody

Catalog Number: 21071-1-AP

1 Publications



## Basic Information

Catalog Number:

21071-1-AP

Size:

600 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15214

GenBank Accession Number:

BC025722

GeneID (NCBI):

136

UNIPROT ID:

P29275

Full Name:

adenosine A2b receptor

Calculated MW:

332 aa, 36 kDa

Observed MW:

36 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse

Cited Species:

mouse

Positive Controls:

WB : mouse brain tissue,

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Hu Guo	39430133	J Orthop Translat	WB,IF

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 21071-1-AP (ADORA2B antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.