

For Research Use Only

ADORA2B Polyclonal antibody, PBS Only

Catalog Number: 21071-1-PBS



Basic Information

| | | |
|---|--|--|
| Catalog Number: 21071-1-PBS | GenBank Accession Number: BC025722 | Purification Method: Antigen affinity purification |
| Size: 1 mg/ml | GeneID (NCBI): 136 | |
| Source: Rabbit | UNIPROT ID: P29275 | |
| Isotype: IgG | Full Name: adenosine A2b receptor | |
| Immunogen Catalog Number: AG15214 | Calculated MW: 332 aa, 36 kDa | |
| | Observed MW: 36 kDa | |

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 21071-1-AP (ADORA2B antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 21071-1-PBS in a different storage buffer formulation.