

For Research Use Only

TTL5 Polyclonal antibody

Catalog Number: 21103-1-AP



Basic Information

Catalog Number:

21103-1-AP

Concentration:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15373

GenBank Accession Number:

BC002766

GeneID (NCBI):

23093

UNIPROT ID:

Q6EMB2

Full Name:

tubulin tyrosine ligase-like family,
member 5

Calculated MW:

1281 aa, 144 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data