

For Research Use Only

# TTL5 Polyclonal antibody, PBS Only

Catalog Number: 21103-1-PBS



## Basic Information

Catalog Number:

21103-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15373

GenBank Accession Number:

BC002766

GeneID (NCBI):

23093

UNIPROT ID:

Q6EMB2

Full Name:

tubulin tyrosine ligase-like family,  
member 5

Calculated MW:

1281 aa, 144 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data