

For Research Use Only

CCDC9 Polyclonal antibody, PBS Only

Catalog Number: 21104-1-PBS



Basic Information

Catalog Number:

21104-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15374

GenBank Accession Number:

BC002787

GeneID (NCBI):

26093

UNIPROT ID:

Q9Y3X0

Full Name:

coiled-coil domain containing 9

Calculated MW:

531 aa, 60 kDa

Observed MW:

60-70 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

CCDC9 is Coiled-coil domain-containing protein 9 which functions on RNA binding. The MW of CCDC9 is about 60-70 kDa in Western blot detection.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 21104-1-AP (CCDC9 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 21104-1-PBS in a different storage buffer formulation.



HeLa cells were subjected to SDS PAGE followed by western blot with 21104-1-AP (CCDC9 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 21104-1-PBS in a different storage buffer formulation.