

For Research Use Only

SMCR8 Polyclonal antibody

Catalog Number: 21125-1-AP

1 Publications



Basic Information

Catalog Number:

21125-1-AP

Size:

640 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15516

GenBank Accession Number:

BC142680

GeneID (NCBI):

140775

UNIPROT ID:

Q8TEV9

Full Name:

Smith-Magenis syndrome
chromosome region, candidate 8

Calculated MW:

937 aa, 105 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
John L Goodier	32678027	Acta Neuropathol Commun	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

