

For Research Use Only

HPSE2 Polyclonal antibody

Catalog Number: 21128-1-AP

1 Publications



Basic Information

Catalog Number:

21128-1-AP

Concentration:

300 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15526

GenBank Accession Number:

BC112356

GeneID (NCBI):

60495

UNIPROT ID:

Q8WWQ2

Full Name:

heparanase 2

Calculated MW:

592 aa, 67 kDa

Observed MW:

54 kDa, 66 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : L02 cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
James Hopkins	39772142	Viruses	WB,IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



L02 cells were subjected to SDS PAGE followed by western blot with 21128-1-AP (HPSE2 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.