

For Research Use Only

HPSE2 Polyclonal antibody

Catalog Number: 21128-1-AP



Basic Information

Catalog Number: 21128-1-AP	GenBank Accession Number: BC112356	Purification Method: Antigen affinity purification
Size: 300 µg/ml	GeneID (NCBI): 60495	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q8WWQ2	
Isotype: IgG	Full Name: heparanase 2	
Immunogen Catalog Number: AG15526	Calculated MW: 592 aa, 67 kDa	
	Observed MW: 54 kDa, 66 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : L02 cells,
Species Specificity: human, mouse, rat	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



L02 cells were subjected to SDS PAGE followed by western blot with 21128-1-AP (HPSE2 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.