

SCGBL Polyclonal antibody

Catalog Number: 21142-1-AP

Basic Information

Catalog Number:

21142-1-AP

Size:

260 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15610

GenBank Accession Number:

BC093909

GeneID (NCBI):

284402

UNIPROT ID:

Q4G0G5

Full Name:

secretoglobin-like

Calculated MW:

96 aa, 11 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IHC 1:20-1:200

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

IHC : human brain tissue, human gliomas tissue, human prostate hyperplasia tissue

Background Information

Storage

Storage:


Store at -20°C. Stable for one year after shipment.

Storage Buffer:


PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human brain using 21142-1-AP (SCGBL antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human brain using 21142-1-AP (SCGBL antibody) at dilution of 1:50 (under 40x lens).