

For Research Use Only

SERPINB7 Polyclonal antibody, PBS Only

Catalog Number: 21170-1-PBS



Basic Information

| | | |
|---|--|--|
| Catalog Number: 21170-1-PBS | GenBank Accession Number: BC069417 | Purification Method: Antigen affinity purification |
| Source: Rabbit | GeneID (NCBI): 8710 | |
| Isotype: IgG | UNIPROT ID: O75635 | |
| Immunogen Catalog Number: AG15394 | Full Name: serpin peptidase inhibitor, clade B (ovalbumin), member 7 | |
| | Calculated MW: 380 aa, 43 kDa | |
| | Observed MW: 43-45 kDa | |

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

Serine protease inhibitor B7 (SERPINB7) is also named Serpin B7, Megsin and TP55. SerpinB7 is identified as a skin-specific endogenous protease inhibitor. SerpinB7 is upregulated in diabetic nephropathy and induces progressive mesangial expansion (PMID: 11877466). SERPINB7 is abundantly expressed in keratinocytes, with particularly high levels in the stratum corneum and stratum granulosum of the epidermis (PMID: 41469525). SERPINB7 promotes keratinocyte differentiation and dampens inflammation in psoriasis in a calcium-dependent manner (PMID: 35864103). SERPINB7 is expressed mostly in renal mesangial cells, and its expression is significantly increased in renal tissues of IgA and diabetic nephropathies. SERPINB7 regulates cell proliferation, signal transduction, and extracellular matrix metabolism (PMID: 37809412).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 21170-1-AP (SERPINB7 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 21170-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 21170-1-AP (SERPINB7 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 21170-1-PBS in a different storage buffer formulation.