

For Research Use Only

ZNF230 Polyclonal antibody

Catalog Number: 21198-1-AP



Basic Information

Catalog Number: 21198-1-AP	GenBank Accession Number: BC050340	Purification Method: Antigen affinity purification
Size: 650 µg/ml	GeneID (NCBI): 7773	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	UNIPROT ID: Q9UIEO	
Isotype: IgG	Full Name: zinc finger protein 230	
Immunogen Catalog Number: AG15551	Calculated MW: 474 aa, 55 kDa	
	Observed MW: 56 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse liver tissue,
Species Specificity: human, mouse	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 21198-1-AP (ZNF230 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.