

For Research Use Only

TXNDC6 Polyclonal antibody

Catalog Number: 21235-1-AP



Basic Information

Catalog Number: 21235-1-AP	GenBank Accession Number: BC107054	Purification Method: Antigen affinity purification
Source: Rabbit	GeneID (NCBI): 347736	Recommended Dilutions: WB: 1:500-1:1000
Isotype: IgG	UNIPROT ID: Q86XW9	
Immunogen Catalog Number: AG15719	Full Name: thioredoxin domain containing 6	
	Calculated MW: 263 aa, 30 kDa	
	Observed MW: 30 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : K-562 cells,
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



K-562 cell lysates were subjected to SDS PAGE followed by western blot with 21235-1-AP (TXNDC6 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.