

For Research Use Only

HMGB3L1 Polyclonal antibody

Catalog Number: 21246-1-AP **1 Publications**



Basic Information

Catalog Number: 21246-1-AP	GenBank Accession Number: BC103851	Purification Method: Antigen affinity purification
Size: 700 µg/ml	GeneID (NCBI): 128872	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: high-mobility group box 3-like 1	
Isotype: IgG	Calculated MW: 130 aa, 15 kDa	
Immunogen Catalog Number: AG15748	Observed MW: 28-31 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, human placenta tissue
Cited Applications: WB	
Species Specificity: human	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xudong Chen	25095979	Tumour Biol	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 21246-1-AP (HMGB3L1 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.