

For Research Use Only

CRYGN Polyclonal antibody

Catalog Number: 21281-1-AP



Basic Information

Catalog Number:

21281-1-AP

Size:

260 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15795

GenBank Accession Number:

BC100880

GeneID (NCBI):

155051

UNIPROT ID:

Q8WXF5

Full Name:

crystallin, gamma N

Calculated MW:

182 aa, 21 kDa

Observed MW:

17 kDa, 25 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse eye tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse eye tissue were subjected to SDS PAGE followed by western blot with 21281-1-AP (CRYGN Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.