

For Research Use Only

SLC39A2 Polyclonal antibody, PBS Only

Catalog Number: 21287-1-PBS



Basic Information

Catalog Number:

21287-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15807

GenBank Accession Number:

BC096723

GeneID (NCBI):

29986

UNIPROT ID:

Q9NP94

Full Name:

solute carrier family 39 (zinc transporter), member 2

Calculated MW:

309 aa, 33 kDa

Observed MW:

33 kDa, 66-70 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Background Information

SLC39A2 (hZIP2) is a human zinc importer in the ZIP family, which is essential in the maintenance of intracellular zinc homeostasis and is upregulated by zinc depletion. SLC39A2 is highly expressed in prostate epithelial cells and is involved in prostate cancer development. Most SLC39 proteins consist of eight transmembrane domains (TMDs), with extracellular or intravesicular amino and carboxy termini (PMID: 34790705). Zip2 gene knockout exacerbated myocardial I/R injury, whereas Zip2 gene transfer reduced myocardial infarction (PMID: 31095941). The predicted molecular weight of SLC39A2 is 33 kDa, and a band of dimers at 66 kDa can also be detected by this antibody (PMID: 24619115).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded Colorectal cancer slide using 21287-1-AP (SLC39A2 antibody) at dilution of 1:100 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 21287-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 21287-1-AP (SLC39A2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 21287-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded Colorectal cancer slide using 21287-1-AP (SLC39A2 antibody) at dilution of 1:100 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 21287-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 21287-1-AP (SLC39A2 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 21287-1-PBS in a different storage buffer formulation.