

For Research Use Only

TMEM117 Polyclonal antibody

Catalog Number: 21314-1-AP

Featured Product

1 Publications



Basic Information

Catalog Number:

21314-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15966

GenBank Accession Number:

BC030284

GeneID (NCBI):

84216

UNIPROT ID:

Q9HOC3

Full Name:

transmembrane protein 117

Calculated MW:

514 aa, 60 kDa

Observed MW:

45-50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : HepG2 cells, MDA-MB-453s cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Meijian Liao	34804929	Front Oncol	WB,IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 21314-1-AP (TMEM117 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.