

For Research Use Only

ZNF23 Polyclonal antibody, PBS Only

Catalog Number: 21322-1-PBS



Basic Information

Catalog Number:

21322-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14755

GenBank Accession Number:

BC007256

GeneID (NCBI):

7571

UNIPROT ID:

P17027

Full Name:

zinc finger protein 23 (KOX 16)

Calculated MW:

643 aa, 73 kDa

Observed MW:

73 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IF/ICC, IP, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 21322-1-AP (ZNF23 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 21322-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of HepG2 cells, using ZNF23 antibody 21322-1-AP at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). This data was developed using the same antibody clone with 21322-1-PBS in a different storage buffer formulation.



IP result of anti-ZNF23 (IP:21322-1-AP, 4ug; Detection:21322-1-AP 1:800) with mouse brain tissue lysate 4000ug. This data was developed using the same antibody clone with 21322-1-PBS in a different storage buffer formulation.