

For Research Use Only

UHRF1 Polyclonal antibody, PBS Only

Catalog Number: 21402-1-PBS

Featured Product



Basic Information

Catalog Number:

21402-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16111

GenBank Accession Number:

BC113875

GeneID (NCBI):

29128

UNIPROT ID:

Q96T88

Full Name:

ubiquitin-like with PHD and ring finger domains 1

Calculated MW:

806 aa, 91 kDa

Observed MW:

91-100 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IF/ICC, IP, Indirect ELISA

Species Specificity:

human

Background Information

UHRF1, also named as ICBP90, NP95 and RNF106, is a putative E3 ubiquitin-protein ligase. It may participate in methylation-dependent transcriptional regulation. UHRF1 binds to inverted 5'-CCAAT-3' box 2 in the TOP2A promoter, and activates TOP2A expression. It is important for G1/S transition and may be involved in DNA repair and chromosomal stability. UHRF1's ability to repress its direct target gene expression is dependent on PHD(UHRF1) binding to unmodified H3R2. In addition, UHRF1 is a novel metabolic guardian restricting AMPK activity(PMID: 34907338).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 21402-1-AP (UHRF1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 21402-1-PBS in a different storage buffer formulation.



WB result of UHRF1 antibody (21402-1-AP; 1:600; incubated at room temperature for 1.5 hours) with sh-Control and sh-UHRF1 transfected HEK-293T cells. This data was developed using the same antibody clone with 21402-1-PBS in a different storage buffer formulation.



IP result of anti-UHRF1 (IP:21402-1-AP, 4ug; Detection:21402-1-AP 1:1000) with HeLa cells lysate 3200ug. This data was developed using the same antibody clone with 21402-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using UHRF1 antibody (21402-1-AP) at dilution of 1:200 and CoralLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 21402-1-PBS in a different storage buffer formulation.