

For Research Use Only

WNT3A Polyclonal antibody

Catalog Number: 21414-1-AP

3 Publications



Basic Information

Catalog Number:

21414-1-AP

Size:

207 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15834

GenBank Accession Number:

BC103921

GeneID (NCBI):

89780

UNIPROT ID:

P56704

Full Name:

wingless-type MMTV integration site family, member 3A

Calculated MW:

352 aa, 39 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human, mouse, rat

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Mengyu Wu	39372681	Stem Cells Int	WB
Shi-Hua Wang	39347028	J Dent Sci	WB
Yiyang Du	37739105	J Ethnopharmacol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

