

IKBKE Polyclonal antibody

Catalog Number: 21431-1-AP

Basic Information

Catalog Number:

21431-1-AP

Size:

213 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16037

GenBank Accession Number:

BC107812

GeneID (NCBI):

9641

UNIPROT ID:

Q14164

Full Name:

inhibitor of kappa light polypeptide
gene enhancer in B-cells, kinase
epsilon

Calculated MW:

716 aa, 80 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

IKBKE is a serine/threonine kinase that plays an essential role in regulating inflammatory responses to viral infection, through the activation of the type I IFN, NF-kappa-B and STAT signaling. It is also involved in TNFA and inflammatory cytokines, like Interleukin-1, signaling.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

