

For Research Use Only

ASB6 Polyclonal antibody

Catalog Number: 21449-1-AP **2 Publications**



Basic Information

Catalog Number:

21449-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13863

GenBank Accession Number:

BC001719

GeneID (NCBI):

140459

UNIPROT ID:

Q9NWX5

Full Name:

ankyrin repeat and SOCS box-containing 6

Calculated MW:

421 aa, 47 kDa

Observed MW:

47 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : L02 cells, HEK-293 cells

Background Information

Notable Publications

| Author | Pubmed ID | Journal | Application |
|-------------|-----------|---------------------|-------------|
| Liyan Gong | 34164402 | Front Cell Dev Biol | WB |
| Qingyong Hu | 38577591 | J Cancer | WB |

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



L02 cells were subjected to SDS PAGE followed by western blot with 21449-1-AP (ASB6 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.